FluoroCore® 2+

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USA: Rx only.

1. PRODUCT DESCRIPTION

FluoroCore® 2+ Dual Cure Core Build-Up Material + Fast Endodontic Post Cement consists of two components, base and catalyst, which when mixed form a dual-cured, highly filled, highly filled, composite resin core build-up and post cementation material. FluoroCore® 2+ Material uses a biocompatible urethane resin and is available in two shades, blue and tooth colored. The blue shade permits clear differentiation of the core material from tooth structure. The tooth colored shade can be used where esthetics and show-through of the core are of primary concern.

1.1 Delivery Forms

FluoroCore® 2+ Material is available in:

- A convenient dual-barelled syringe. This delivery system simplifies direct, precise intraoral delivery and minimizes product waste.
- Two shades, blue and tooth colored

1.2 Composition

Urethane Dimethacrylate
Di- & Tri-functional Methacrylates
Barium Boron Fluoroaluminosilicate Glass
Camphorquinone (CQ) Photoinitiator
Photoaccelerators
Silicon Dioxide
Benzoyl Peroxide

1.3 Indications

- Vital or non-vital tooth core build-up (replacement of existing restorations and/or lost tooth structure) as a base prior to fabricating an indirect restoration.
- Cementation of endodontic posts

1.4 Contraindication

FluoroCore® 2+ Material is contraindicated in patients with a history of allergic reaction to methacrylate resins or any other of the components.

1.5 Compatible adhesives

FluoroCore® 2+ Material is to be used after the application of a suitable dentin and enamel adhesive. FluoroCore® 2+ Material is compatible with numerous dentin/enamel adhesive systems including all DENTSPLY adhesives designed for use with dual/self-cured composite restoratives (see complete Directions for Use of selected adhesive). Other dentin/enamel adhesive systems may be used at the discretion and sole responsibility of the dental practitioner.

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these directions for use.

2.1 Warnings

FluoroCore® 2+ Material contains polymerizable methacrylates which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.

- Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.

2.2 Precautions

- This product is intended to be used only as specifically outlined in the Directions For Use. Any use of this product inconsistent with the Directions For Use is at the discretion and sole responsibility of the practitioner.
- Contact with saliva and blood during composite placement may cause failure of the restoration. Use of rubber dam or adequate isolation is recommended.
- Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
- Syringe should be tightly closed by replacing the original cap immediately after use.
- FluoroCore® 2+ Material is intended to be used as a dual-cure product. To obtain optimal curing follow the Step-by-Step instructions. Self-cured, dual-cured and light-cured methods are integrated into the procedural steps.
- FluoroCore® 2+ Material should be mixed in equal volumes directly through the supplied mix tips for optimum performance. Variations from this 1:1 ratio may adversely affect working time and compromise final physical properties. Dispensing and mixing other than through the supplied mixing tips is not recommended.
- FluoroCore® 2+ Material should extrude easily. DO NOT USE EXCESSIVE FORCE. Excessive pressure may result in unanticipated extrusion of the material or cause syringe rupture.
- FluoroCore® 2+ Material when used for post cementation has a working time of *only 40 seconds!* Prepare and organize the necessary materials and instruments to allow a fast insertion of the post.
- Interactions
 - Eugenol-containing dental materials should not be used with FluoroCore® 2+ Material because they may interfere with hardening and cause softening of the polymeric components of the materials.
 - Variable in-vitro data exist regarding use of light-cured-only adhesives without Self Cure Activator in conjunction with dual-cured or self-cured composite restorative materials such as FluoroCore® 2+ Material in limited or no light curing applications. Chemical/Product incompatibility may adversely affect product efficacy, leading to premature restoration failure.
 - As with any dual cure resin system, the use of a dual cure adhesive system can shorten working time. This effect should be investigated in the laboratory prior to clinical use.
 - FluoroCore® 2+ Material is a light-cured material. Therefore, it should be protected from ambient light. Proceed immediately once material has been placed.

2.3 Storage conditions

FluoroCore® 2+ Material syringes should be stored with the original cap only, at refrigerated temperatures between 2°C/35°F and 8°C/46°F. Keep out of direct sunlight and protect from moisture. Allow material to reach room temperature prior to use. Do not freeze. Do not use after expiration date. Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

2.4 Adverse reactions

Product may irritate the eyes and skin. Eye contact: irritation and possible corneal damage. Skin contact: irritation or possible allergic response. Reddish rashes may be seen on the skin. Mucous Membranes: inflammation (See Warnings).

3. STEP-BY-STEP INSTRUCTIONS

Working and Setting Times for Core Build-up and Endodontic Post Cementation

Maximum Working Time at body temperature (37°C) 40 seconds

for manipulation/post seating after intraoral placement

Minimum Working Time at room temperature (22°C) 1 minute 30 seconds

for application after start of mix

Setting Time Light Cure 20 seconds' (3mm depth)

Setting Time Self Cure (catalyst initiated) (37°C) 3 minutes

3.1 Core build-up

3.1.1 Cleaning and Preparation of Tooth

- Finish preparation by removing existing restorations and/or caries.
- Clean freshly cut enamel and dentin with water spray and air dry. Do not desiccate dentin
- Place pin, post or matrix as needed.

3.1.2 Pulp protection, Tooth Conditioning/Dentin Pretreatment, Adhesive Application
Refer to applicable manufacturer's directions for pulp protection, tooth conditioning and/or adhesive application. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of FluoroCore® 2+ Material.

3.1.3 Dispensing and mixing FluoroCore® 2+ Dual barreled syringe

- 

CAUTION

Danger of injury due to excessive force

 - Apply slow and steady pressure on the syringe plunger
 - Do not use excessive force — syringe rupture may result

- Remove from refrigerator and bring to room temperature
- Remove syringe cap and dispense and discard a small amount of material from the dual barreled syringe. Be sure material is flowing freely from both ports. Holding syringe vertically, carefully wipe away excess so base and catalyst do not cross contaminate and cause obstruction of the ports. Save syringe cap for replacement following use.
- Install a mixing tip, with optional intraoral tip attached if desired, on the cartridge by lining up the v-shaped notch on the outside of the mixing tip with the v-shape notch on the syringe flange. Turn colored mixing tip cap 90 degrees in a clockwise direction to lock in place on syringe.
- Gently depress syringe plungers to begin the flow of material. **DO NOT USE EXCESSIVE FORCE.** If force is encountered, remove syringe from operating field, remove and discard mixing tip. Check for obstruction and confirm material is flowing from both syringe barrels. Wipe barrels and install new mixing tip as outlined above. Dispense a small amount through the mixing tip onto a mixing pad and discard.
- Proceed immediately to placement, Step 3.1.4. Placement of FluoroCore® 2+ Material.

3.1.4 Placement of FluoroCore® 2+ Material

- FluoroCore® 2+ Material produces a syringable mixed consistency for direct placement using the dual-barreled syringe and a banded matrix retainer or clear crown form. Alternatively, mixed material may be expressed onto a clean paper pad and manually placed with the operator’s placement instrument(s) of choice. Working time of mixed material at room temperature (22°C/72°F) is approximately 1’30” (one minute thirty seconds) when protected from ambient light. Elevated oral temperatures and/or contact with dual cure adhesives may shorten work time (see 2.2 Precautions). Any of the following placement techniques are acceptable. The preferred method is to use both visible light and chemical self cure (see Step 3.1.4.3 Dual Cure).
- Visible Light Cure:* Place mixed material directly in 2-3mm increments and light cure increment for 20 seconds’.
 - Chemical Self-Cure:* Place large increment(s) or load mixed material into a clear crown form and seat onto the preparation. Allow to self-cure approximately 3 minutes. Visible light curing’ of external surfaces is optional, but advisable.
 - Dual Cure:* After placing the first 2-3mm increment, and visible light curing for 20 seconds’, large increments may then be placed or material loaded into a clear crown form, seated onto the preparation and allowed to self cure. Visible light curing of external surface is optional, but advisable.

3.1.5. Finishing the FluoroCore® 2+ Material

After initial self cure set, and removal of matrix, a surface light-cure of the material for an additional 20 seconds’ is recommended. Initial crown preparation (gross reduction) may begin immediately using rotary instruments. If used per Step 3.1.4.2 (Chemical Self-Cure) or 3.1.4.3 (Dual Cure), utilizing chemical self cure, greatest bulk hardness is obtained 6 minutes after start of mixing at which time the final preparation may be completed.

3.2 Endodontic Post Cementation

Tooth Conditioner/Dentin Pretreatment/Adhesive application

Refer to applicable manufacturer’s directions for tooth conditioning and/or adhesive application. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of FluoroCore® 2+ Material, as outlined below. **NOTE:** If core build-up will immediately follow post cementation, treat all available tooth structure beyond post space with adhesive for core buildup bonding.

- Preparation of post: Treat surface of post according to adhesive manufacturer’s and/or dental laboratory’s instructions, i.e., etching, mechanical roughening, adhesive application.

- Dispensing and mixing FluoroCore® 2+ Dual barreled syringe: Prepare syringe, bleed material and install mixing tip (with endo intraoral tip attached) as outlined above in 3.1.3. (Dispensing and mixing FluoroCore® 2+ dual-barreled syringe, Steps 1-4)

- Proceed immediately with placement. Dispense the same material directly from the syringe onto the post surface and into the post space through the orifice of the post hole preparation.

- Immediately place post into final position and stabilize. **IMPORTANT NOTE:** Work time of applied cement within a post space pre-treated with dual cure adhesive is 40 seconds maximum from application. Post must be fully seated and stabilized within 40” (see 2.2 Precautions, Interactions sections).

- 

CAUTION

Inadequate insertion of post due to short working time

 - At body temperature working time is 40 seconds from start of application
 - Prepare and organize necessary material and instruments to allow a fast insertion of the post

Once stabilized, light cure, using a visible light, curing unit¹. Cure all accessible areas of the post for 20 seconds.

- 

CAUTION

Inadequate polymerization due to insufficient curing

 - Check compatibility of curing light
 - Check curing cycle
 - Check curing output before each procedure

- Proceed with the core build-up as described in Section 3.1. (Core Build-Up)

4. HYGIENE

4.1 Cleaning the syringe

For dual-barreled syringe, remove used mixing tip and discard appropriately. Replace original syringe cap prior to storing. The dual-barreled syringe may be cleaned by scrubbing with a disposable towel soaked with hot water and soap or detergent.

4.2 Disinfecting the syringe

Disinfect the syringe with a hospital-level, tuberculocidal disinfectant solution according to national/local regulations. Iodophors, sodium hypochlorite (5.25%), chlorine dioxide and dual or synergized quaternary ammoniums are approved disinfectants. Some phenolic-based agents and iodophor-based products may cause surface staining. The disinfectant manufacturer’s directions should be followed properly for optimum results. Water-based disinfectant solutions are preferred.

5. LOT NUMBER AND EXPIRATION DATE

- Do not use after expiration date. ISO standard uses: “YYYY/MM.”
- The following numbers should be quoted in all correspondence:
 - Reorder Number
 - Lot number
 - Expiration Date

^[1] Use suitable visible light curing unit designed to cure materials containing camphorquinone (CQ) initiator, i.e. spectral output containing 470 nm. Minimum light output must be at least 550mW/cm² exposure. Refer to curing light manufacturer’s recommendations for compatibility and curing recommendations.

©2009 DENTSPLY International. All Rights Reserved. Printed in U.S.A.

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DENTSPLY DeTrey GmbH